


TECHNICAL SHEET

Reference	Oral tip syringe			
Trade Name:	Oral tip syringe			
Product Description	Oral tip syringes are used to deliver medication and feed orally. Used by patients in the hospital and community environment.			
Intended Use	To deliver feed/Medication orally			
Raw Materials:	<ul style="list-style-type: none"> • COMPONENT 1: Barrel-Polypropylene • COMPONENT 2: Plunger-Polypropylene • COMPONENT 3: Sealing ring-Silicon rubber • COMPONENT 4: Lubricant-Medical grade silicone oil • COMPONENT 5: Oral tip - Polypropylene 			
				
Packaging	Multiple package (see table below). QML (minimum delivery quantity): 1 Box			
	CODICE	ML	PACKAGE Pcs.	BOX Pcs.
	OT005	0.5ml	1	100
	OT01	1ml	1	100
	OT05	5ml	1	100
	OT10	10ml	1	100
	OT20	20ml	1	80
	OT25	2.5ml	1	100
Shelf Life:	5 years from the date of sterilization			
Usage	Single use device			
Storage	Normal storage procedures and conditions			
Disposal	The user must follow the legal regulations regarding disposal of hospital waste.			
Manufacturer	Medicina ltd Units 1-4 Rivington View Business park Station Road Blackrod Bolton BL6 5BN			
Classification	Class IIa Annex V & VII			
Conformity Assessment Route:	MDD 93/42/EEC			