




OT01LD	1ml low dose oral syringe
OT25LD	2.5ml low dose oral syringe

Medicina low dose oral syringe (ENFit compatible)

Medicina Ltd. has developed a low dose syringe that is designed to give medication or enteral feed orally.

The oral tip one piece syringe has a smooth barrel and tip to make giving oral medication and feeds less traumatic for small babies and children.

Cup fill is not recommended. If used please remove any fluid from the "moat" this is not in the fluid path but will go into the threads of the male ENFit connector and may cause connection problems



Oral tip syringes




OT005 0.5ml oral tip syringe (100)




OT01 1ml oral tip syringe (100)




OT01LD 1ml low dose oral tip syringe* (100)




OT25 2.5ml oral tip syringe (100)




OT25LD 2.5ml low dose oral tip syringe* (100)



OT05 5ml oral tip syringe (100)



OT10 10ml oral tip syringe (100)



OT20 20ml oral tip syringe (100)

Syringe caps




OTC01 Oral tip syringe cap (50)

*Please note OT01LD and OT25LD low dose syringes require the M158LD caps


ENFit syringes




LPE01LD 1ml low dose enteral syringe (100)



LPE25LD 2.5ml low dose enteral syringe (100)




LPE05 5ml enteral syringe (100)




LPE10 10ml enteral syringe (100)




LPE20 20ml enteral syringe (80)



LPE30 30ml enteral syringe (80)



LPE60 60ml enteral syringe (60)



LPE100 100ml enteral syringe (30)



M158/1LD Male Caps (single) (50)

M158/8LD Male Caps
50 packs of 8 (400)

Syringe accessories



OTC01 Oral tip syringe cap (50)



BA01L Stepped Bottle Adapter (50)



BA02L 15-16.5mm Bottle Adapter (50)

BA03L 18-19.5mm Bottle Adapter (50)

BA04L 25.5-27mm Bottle Adapter (50)

BA05L 18.5-20.5mm Bottle Adapter (50)
(Oramorph)

BA06L 20-21.5mm Bottle Adapter (50)

BA07L 10.5-11mm Bottle Adapter (50)

BA08L 15.5-18mm Bottle Adapter (50)



RC075L 75mm Rectal Medicine
Straw (50)

RC150L 150mm Rectal Medicine
Straw (50)



C005L 5cm Medicine Straw (50)



C020L 20cm Milk Straw (50)

C014L 14cm Milk Straw (50)

C010L 10cm Milk Straw (50)



MF01L Filter Straw (50)



at the heart of healthcare

