



SURGICAL GLOVE SOLUTIONS DEVELOPED, PRODUCED & PACKAGED IN EUROPE



MILE STONES

GROUNDBREAKING BEGINNINGS

Semperit is founded by Johann Nepomuk Reithoffer, a tailor, entrepreneur and inventor specialising in waterproof fabrics. His patented weaving machine marks the start of the European rubber industry.

OUR FUTURE HOME

Construction commences on a factory in Wimpassing, Austria, which will one day house Sempermed's surgical glove production and research and development centre.



FIRST GLOVES Semperit starts producing latex gloves.

GLOBAL APPEAL

Sempermed surgical gloves are donned by professionals around the world.

A NEW EUROPEAN MARKET LEADER

Sempermed[®] Supreme, one of the most frequently used powder-free surgical gloves in Europe, establishes itself as a glove of choice for surgeons.

INNOVATION AWARD

Our patented photochemical process, used to manufacture the skin-friendly, accelerator-free Sempermed® Syntegra UV, is awarded the prestigious EARTO innovation award.

SUSTAINABILITY & SKIN-FRIENDLINESS

The innovative Green Glove Technology introduces a new level of sustainability in glove manufacturing.



A CENTURY OF GLOVE EXPERTISE

Semperit Investments Asia Pte. Ltd, 8 Jurong Town Hall Road, #29–03 to 06 The JTC Summit, Singapore 609434, Singapore Tel.: +65 6274 4861 · Fax.: +65 6274 6977 sempermed@semperitgroup.com · www.sempermed.com Contact Us

Sempermed Europe GmbH Am Belvedere 10, 1100 Vienna, Austria Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630 sempermed@semperitgroup.com · www.sempermed.com

WHY SURGEONS TRUST SEMPERMED

Listening to healthcare professionals and observing them at work, we understand which innovations make a real difference in the operating room. We continuously push for new products that empower and protect surgical staff, care for even the most sensitive skin, and are manufactured in a more sustainable manner.



DOUBLE PROTECTION

Sempermed surgical gloves are double-certified*, in accordance with Medical Device Directive 93/42/EEC as well as Personal Protective Equipment Regulation (EU) 2016/425, to ensure your safety. Thanks to an optimised wall thickness, our gloves offer superior tactile sensitivity without making any concessions when it comes to your protection. *except Sempermed Senso

SAFE FOR SENSITIVE SKIN

We offer solutions that protect your skin from common allergies. Synthetic materials which are free from allergenic latex proteins, skin-friendly manufacturing methods like our patented UV technology, and our specially developed washing and finishing process significantly reduce the risk of skin allergies.

QUICK AND EASY TO DON Quick don

Our surgical gloves possess excellent donning properties thanks to the specially developed Quick Don[™] technology. The powder-free synthetic inner coating reduces friction resistance, making our easyto-don surgical gloves ideal for double-donning and intraoperative glove changes.

A COMPREHENSIVE PORTFOLIO

Surgeons and OR personnel encounter highly sensitive or high-risk procedures in their daily work. Our extensive portfolio, containing features like a micro-rough outer surface or an optimised wall thickness, offer solutions that empower and protect surgical staff and patients.

FULLY ANATOMICAL SHAPE

With a design based on a natural, relaxed hand posture featuring a slightly wider back of the hand, our gloves allow for superior mobility and ensure that hands experience significantly less fatigue, especially during long procedures.

MADE IN AUSTRIA SINCE 1920

We develop and produce surgical gloves in Austria. As a globally operating company with strong European values, we are committed to the International Labour Organisation's standards, are an amfori Business Social Compliance Initiative member, and are always on the lookout for more sustainable solutions.

IMPORTANT NOTE: The latest produ Products' tab. Failure to observe this frequency of use and tolerability of th Semperit is not liable for incorrect storage and handling. In case of da used. In case of any doubt, expert ad cation is up to date as of the date of p and may vary in individual cases. This d CAUTION: Natural latex can cause all no liability for allergic reactions. "Co ation is available at www.sempermed.com under the on, in particular with regard to (chemical) resistance, can result in personal injury and/or material damage. e gloves nor will Semperit be liable for improper unusual signs of wear and tear, gloves shall not be d be sought before use. The information and classifi-. Technical details are average values from production s subject to mistakes, printing errors and amendments. tions, including anaphylactic shock. Semperit bears 21 Sempermed

	NATURAL LATEX FREE					
BRAND NAME	sempermed°() syntegra UV	semper <mark>med</mark> °() syntegra IR	semper med ®() syntegra greer			
	ACCEL RECCEL RECCEL RECCEL RECERTOR	NATURAL REE Cuick don **				
PRODUCT BENEFITS	AWARD-WINNING PROTECTION FOR YOUR SKIN	ADVANCED SKIN PROTECTION & PERFORMANCE	THE INNER GLOVE FOR HIGH- RISK PROCEDURES			
	 Polyisoprene: excellent comfort & protection Low-allergenic ¹: accelerator-free & latex-free Great comfort: extraordinary latex-like fit & feel Fast & easy to don: Quick Don[™] inner coating 	 Reduced allergy risk: innovative accelerators, latex-free Fast & easy to don: Quick Don™ inner coating Great comfort: extraordinary latex-like fit & feel Anatomical shape: reduces strain for wearer 	 Ideal inner glove: designed for double-donning Skin-friendly: innovative accelerators, latex-free Great comfort: extraordinary latex-like fit & feel Wall thickness optimal for high-risk surgeries 			
ТҮРЕ	Sterile surgical glove for single use, powder-free, with synthetic lining					
MATERIAL	Synthetic polyisoprene	Synthetic polyisoprene	Synthetic polyisoprene			
COLOUR	White	Cream	Green			
GLOVE SHAPE	Fully anatomical with curved fingers and rolled rim					
SIZE/OVERALL LENGTH as per EN 455-2	5 ½, 6 and 6½ median 270mm 7, 7½ and 8 median 280mm 8½ and 9 median 285mm					
WALL THICKNESS measured in single layer in the palm area	0.20-0.25 mm	0.19-0.24 mm	0.19–0.24 mm			
BARRIER PERFORMANCE	AQL 0.65					
FORCE AT BREAK as per EN 455-2	≥9N					
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	3 years					
STERILISATION		Radiation STERILE R				
PACKAGING	Left and right glove turned up cuff in inner pouch, ozone-tight, sealed in flat medical peel pack. In dispenser carton with sterilisation indicator: 40 pairs. In transport carton with sterilisation indicator: 240 pairs.					
PURPOSE according to MP Regulation 93/42/EWG / according to PPE Regulation (EU) 2016/425		Medical device class IIa ingle-use protective glove – PPE category ion against chemical substances for limited				
ARTICLE NUMBERS	EU Size 5 ½ 827058521 Size 6 827058601 Size 6½ 827058621 Size 7 827058701 Size 7½ 827058721 Size 8 827058801 Size 8½ 827058821	EU Size 5 ½ 827056521 Size 6 827056601 Size 6 ½ 827056621 Size 7 827056701 Size 7 ½ 827056721 Size 8 827056801 Size 8 ½ 827056821	EU Size 5 ½ 827053521 Size 6 827053601 Size 6 ½ 827053621 Size 7 827053701 Size 7 ½ 827053721 Size 8 827053801 Size 8 ½ 827053821			
	Size 9 827058901	Size 9 827056901	Size 9 827053901			

¹⁾ Crepy et.al. (2018), Accelerator-free gloves as alternatives in cases of glove allergy in healthcare workers, in Contact Dermatitis. 2018 Jan;78(1):28-32. doi: 10.1111/cod.12860. Epub 2017 Jul 27. https://www.ncbi.nlm.nih.gov/pubmed/28748553 ² Please note the Sempermed Chemical Resistance list, available on www.sempermed.com/sempermed_informs/chemical_resistance

	NATURAL RUBBER LATEX (NR)					
BRAND NAME	semper med®(§ supreme	sempermed® supreme plus	sempermed [®] (s) supreme green			
	quick don	[*] quick d	on [™] quick don [™]			
PRODUCT BENEFITS	A PREMIUM ALL-ROUNDER YOU CAN TRUST	OPTIMAL GRIP IN EVERY SITUATION	DOUBLE PROTECTION THAT YOU CAN TRUST			
	 Extra-supple: great comfort & tactile sensitivity Fast & easy to don: Quick Don[™] inner coating Anatomical shape: reduces strain for wearer Skin-friendly: washed to minimise skin irritants 	 Micro-rough surface: for dry & damp applications Fast & easy to don: Quick Don[™] inner coating Anatomical shape: reduces strain for wearer Skin-friendly: washed to minimise skin irritants 	 Clear detection of perforations: when double-donning with a lighter colored outer glove Tactile sensitivity thanks to optimised wall thickness Ideal inner glove: designed for double-donning As underglove for Sempermed[®] supreme or supreme plus 			
ТҮРЕ	supreme or supreme plus Sterile surgical glove for single use, powder-free, with synthetic lining					
MATERIAL	Natural rubber latex (NR)	Natural rubber latex (NR)	Natural rubber latex (NR)			
COLOUR	Natural white	Natural white	Green			
GLOVE SHAPE	Fully anatomical with curved fingers and rolled rim					
SIZE/OVERALL LENGTH as per EN 455-2		5½, 6 and 6½ median 270mm 7, 7½ and 8 median 280mm 8½ and 9 median 285mm				
WALL THICKNESS measured in single layer in the palm area	0.19–0.23 mm	0.19-0.23 mm	0.17–0.20mm			
BARRIER PERFORMANCE	AQL 0.65					
FORCE AT BREAK as per EN 455-2	\geq 12N, exceeding EN 455-2 (\geq 12N instead of 9N)					
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	3 years					
STERILISATION		Radiation STERILE R				
PACKAGING	Left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in flat medical peel pack. In dispenser carton with sterilisation indicator: 50 pairs. In transport carton with sterilisation indicator: 300 pairs.					
PURPOSE according to MP Regulation 93/42/EWG / according		Medical device class IIa Single-use protective glove – PPE cate	gory III			

Regulation 93/42/EWG / according to PPE Regulation (EU) 2016/425		Single-use protective glove – PPE category III (protection against chemical substances for limited time) ²⁾					
		EU	INT		EU		EU
	Size 51/2	822751521	822751525	Size 51/2	822851521	Size 51/2	823756521
	Size 6	822751601	822751605	Size 6	822851601	Size 6	823756601
	Size 61/2	822751621	822751625	Size 61/2	822851621	Size 61/2	823756621
ARTICLE	Size 7	822751701	822751705	Size 7	822851701	Size 7	823756701
NUMBERS	Size 7 1/2	822751721	822751725	Size 7 1/2	822851721	Size 7 1/2	823756721
	Size 8	822751801	822751805	Size 8	822851801	Size 8	823756801
	Size 81/2	822751821	822751825	Size 8 1/2	822851821	Size 81/2	823756821
	Size 9	822751901	822751905	Size 9	822851901	Size 9	823756901

NATURAL RUBBER LATEX (NR)

semper**med**®(5) supreme duo

semper**med**®(§) derma PF





quick don[™]



BRAND NAME

THE OPTIMAL DOUBLE-DONNING SOLUTION

- Faster to open: two gloves in one peel pack
- (Sempermed[®] supreme + Sempermed[®] supreme green) **Easy to double-don** thanks to an underglove that's a ½ size bigger
- Lasy to double-don thanks to an undergrove that's a 72 size
- Optimal safety: clear detection of perforations
- Wall thickness optimised for great fit & feel



A HIGHLY VALUED PROCEDURE GLOVE

- A premium fit & feel valued by frequent users
- Great value: a quality procedure glove
- Confident grip: due to the micro-rough surface
- Fast & easy to don: Quick Don[™] inner coating

ТҮРЕ	Sterile surgical glove for single use, powder-free, with synthetic lining				
MATERIAL	Natural rubber latex (NR)	Natural rubber latex (NR)			
COLOUR	Natural white + Green	Natural white			
GLOVE SHAPE	Fully anatomical with curved fingers and rolled rim				
SIZE/OVERALL LENGTH as per EN 455-2	5½, 6 and 6½ median 270mm 7, 7½ and 8 median 280mm 8½ and 9 median 285mm	5 ½, 6 and 6 ½ median 260 mm 7, 7 ½ and 8 median 270 mm 8 ½ and 9 median 280 mm			
WALL THICKNESS measured in single layer in the palm area	0.19–0.23mm + 0.17–0.20mm	0.16–0.21 mm			
BARRIER PERFORMANCE	AQL 0.65				
FORCE AT BREAK as per EN 455-2	≥ 12N, exceeding EN 455-2 (≥ 12N instead of ≥ 9N)				
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	3 years				
STERILISATION	Radiation STERILE R				
PACKAGING	Left and right glove turned inside out, in inner pouch, ozone-tight, sealed in medical peel pack, double-packed. In dispenser carton with sterilisation indicator: 25 pouches. In transport carton with sterilisation indicator: 150 pouches.				
PURPOSE according to MP Regulation 93/42/EWG / according to PPE Regulation (EU) 2016/425	Medical device class IIa Single-use protective glove – PPE category III (protection against chemical substances for limited time) ²⁾				
	EU	EU			
ARTICLE NUMBERS	Size 5 ½ 3000007261 Size 6 300007262 Size 6½ 300007263 Size 7 300007264 Size 7½ 300007265 Size 8 300007266 Size 8½ 300007267	Size 5 ½ 822751522 Size 6 822751602 Size 6 ½ 822751622 Size 7 822751702 Size 7 ½ 822751722 Size 8 822751802 Size 8 ½ 822751802 Size 8 ½ 822751802			
	Size 9 3000007268	Size 9 822751902			