















## Wheelchair weigher Conformity assigned

6708

- Conformity assessed according to Class III, MPG Class 1 m
- Platform with integrated drive-on ramps
- Including foldable railing, white painted
- ▶ Transportable due to two integrated castors and handle
- Battery and mains operation







## **Features**

- Conformity assessed according to Class III, MPG Class 1 m
- Platform with integrated drive-on ramps
- Including foldable railing, white painted
- Passage width: 760 mm
- Terminal with plastic housing
- Backlit LC display, digit height: 20 mm
- Transportable due to two integrated castors and handle
- Functions: weighing, taring, hold, BMI
- ▶ Battery level indicator (battery symbol in the display)
- ▶ Energy-saving auto-off function
- Units in kg
- ▶ Battery and mains operation 100 240 V, 50-60 Hz (included in scope of delivery)

## **Products**

Order number	Approvable	<b>Division</b>	<b>Maximum load</b>	<b>Platform size</b>	<b>Product weight</b>
Variation		g	kg	mm	kg
6708 Wheelchair weigher	Yes	100	300	970x840x82	46